



FEATURES

- High resolution test images • ETDRS chart • Single row and simple optotype display of charts
- Fixation targets for children & educational charts • Compact and user friendly remote • Free from interference of other remotes
- Power saving mode • Software upgradable

SPECIFICATIONS



The image shows the AuroChart digital vision chart device. It consists of a monitor displaying the 'auroChart' logo and a keyboard. The device is shown in a perspective view, with a red arrow pointing towards the specifications table.

Display size	: 18.5 inch
Resolution	: 1280 x 720
Luminance	: 250 ± 30 cd/m ²
Power Supply	: 100 - 240 VAC ± 10%, 50-60 Hz UPS connection is recommended
Power Consumption	: 23 Watts
Dimension (mm)	: 449 x 361 x 183 (W x H x D)
Weight	: 2.5 kg
Power Saving Mode	: Activates when no key is pressed for 2 hours
Remote Control	: IR
Mounting	: Wall Mount
Temperature	: 10 - 45 ° C
Humidity	: 10 - 80 % Rh
Testing Distance	: 10 feet (3 meters) and 8 feet (2.4 meters)
Standard Charts	: Snellen, logMAR, ETDRS, Dot acuity test and Paediatric
Special Charts	: Astigmatism, Duo chrome, Amsler grid, Contrast Sensitivity tests, Color Vision Test, Phoria tests, Aniseikonia test, Fixation Targets
Educational Charts	: Structures of the Eye, Refractive errors and Charts for ocular conditions
Languages	: Tamil, Telugu, Kannada, Malayalam, Hindi, Gujarati, Oriya and Punjabi
Accessories	: Remote Control User Manual Installation Kit

* Information published in this catalogue is subject to change without notification



No 1, Sivagangai Main Road, Veerapanjan, Madurai 625 020, India

Phone: 91 452 2446100 . Fax: 91 452 2446200 . E-mail: sales@aurolab.com . Web: www.aurolab.com

