

auro@HART

Digital vision chart



FEATURES

- High resolution test images • ETDRS chart • Single row and simple optotype display of charts
- Fixation targets for children & educational charts • Power saving mode
- Compact and user friendly remote with menu driven options

 **aurolab**
Excellence . . . in sight



SPECIFICATIONS

Display size	: 18.5 inch
Resolution	: 1280 x 720
Luminance	: 200 ± 30 cd/m ²
Power Supply	: 100 - 240 VAC ± 10%, 50-60 Hz UPS connection is recommended
Power Consumption	: 23 W
Dimension (mm)	: 449 x 361 x 183 (W x H x D)
Weight	: 2.5 kg
Remote Control	: IR
Mounting	: Wall Mount
Temperature	: 10 - 45 ° C
Humidity	: 10 - 80 % Rh
Testing Distances	: 2.4 m, 3.0 m, 4.0 m, 5.0 m, 6.0 m & 6.0 m Mirror (Distance selection can be made by accessing menu via remote)
Standard Charts	: Snellen, logMAR, ETDRS, Dot acuity test and Paediatric
Special Charts	: Astigmatism, Duo chrome, Amsler grid, Contrast Sensitivity tests, Color Vision Test, Phoria tests, Aniseikonia test, Fixation Targets
Educational Charts	: Structures of the Eye, Refractive errors and Charts for ocular conditions
Languages	: Tamil, Telugu, Kannada, Malayalam, Hindi, Gujarati, Oriya, Punjabi, Bengali & Urdu (5 languages can be selected at a time by accessing menu via remote)
Accessories	: Remote Control User Manual Key Cards



Information published in this catalogue is subject to change without notification



No 1, Sivagangai Main Road, Veerapanjan, Madurai 625 020, India

Phone: 91 452 2446100 . Fax: 91 452 2446200 . E-mail: info@aurolab.com . Web: www.aurolab.com

