

auroSling

ptosis sling
For perfect ptosis surgery



- **Easy to Bury**
 - 7mm sleeve
- **Convenient**
 - Malleable needles

Specification

Aurosling is an implantable flexible silicone rod attached to malleable sharp needles with a silicone sleeve to provide,

- Excellent biocompatibility
- Good lid closure
- For easy adjustment

Product Code	Needle		Length of Silicone rod	Length of Silicone sleeve
	Length	Diameter		
T9052	6.3cm	920micron	40cm	7mm



 **aurolab**
Excellence . . . in sight

















ARAVIND EYE CARE SYSTEM

ORBIT & OCULOPLASTY SUTURE NEEDLES

aurolac Lacrimal Intubation set

Specification: Lacrimal intubation set comprises of two flexible stainless steel probes attached through a hollow medical graded silicone tube.

Product code	Probe		Silicone tube length - 30 cm	
	Length	Diameter	Silicone tube ID	Silicone tube OD
L6032*	11 cm	0.60 mm (23G)	0.30 mm	0.64 mm

Needle	Suture Material	Suture Length	Suture Sizes					
			9-0	8-0	6-0	5-0	4-0	2-0
Advanced Side Cutting Spatula Diameter - 0.28mm Length - 8.0mm 3/8 Circle - 135° 	ausilk Silk Black Braided Silicone Coated	76cm (30") 			B6105			
Conventional Cutting Diameter - 0.43mm Length - 15.9mm 3/8 Circle - 135° 	ausilk Silk Black Braided Silicone Coated	76cm (30") 					B4401	
Taper Point Diameter - 0.66mm Length - 26.2mm 1/2 Circle - 170° 	aubond Polyester Green Braided	90cm (35") 						E2202
Taper Point Diameter - 0.71mm Length - 17.5mm 1/2 Circle - 180° 	aubond Polyester Green Braided	90cm (35") 						E2204
Reverse Cutting Diameter - 0.38mm Length - 13.1mm 3/8 Circle - 130° 	aurolene Polypropylene Blue Monofilament	70cm (28") 				P5301		
Taper Point Diameter - 0.10mm Length - 5.1mm 3/8 Circle - 135° 	aurolene Polypropylene Blue Monofilament	15cm (6") 	P9201					
Advanced side cutting spatula Diameter - 0.33 mm Length - 8.0 mm 1/4 Circle - 90° 	Polycryl PolyGlycolic Acid Violet Braided	45cm (18") 			A6108*			



Except*

Information published in this catalog is subject to change without notification

 **aurolab**
Excellence ... in sight

For Further Details Contact
Product Manager
Suture - Needle Division
Mobile: 98947 97163
E-mail: v.saravanam_oorthy@aurolab.com

No 1, Sivagangai Main Road, Veerapanjan, Madurai - 625 020, India
Phone: 91 452 2446100 Fax: 91 452 2446200
E-mail: sales@aurolab.com Web: www.aurolab.com

Issue 3 - 08/10



ARAVIND EYE CARE SYSTEM