

# AUROSPHERE

For enucleation & evisceration procedure

Aurosphere is made up of Permanent implantable grade silicone material for excellent biocompatibility.

Aurosphere does not allow direct or indirect integration with the prosthesis. It provides

- Minimal risk of implant migration, exposure or infection
- Inexpensive enucleation option
- Improves both cosmesis and motility.
- Base for prosthesis
- Safe and Reliable

Aurosphere is indicated in all cases of enucleation or evisceration.

## Specification:

Product Code	AS - 14	AS - 16	AS - 18	AS - 20
Diameter	14 mm	16 mm	18 mm	20 mm



# AUROSTENT

For canalicular laceration


Aurostent, a medical graded silicone implant is an effective tool in reconstructing traumatic canalicular lacerations.

- Carries minimum risk of injury to the non-lacerated canaliculus
- Results in high anatomical and functional success rates
- Safe and Cost effective
- Securely anchored at the punctum






















It is an ideal stent which is inert, soft, and pliable to minimize ocular irritation; stable enough to remain in place long enough for pericanalicular inflammatory reactions to subside.

## Specification:

Silicone rod length	40 mm
Silicone rod Diameter	0.64 mm



# SUTURE RANGES FOR ORBIT AND OCULOPLAST

Needle	Suture Material	Suture Length	Suture Sizes						
			9-0	8-0	7-0	6-0	5-0	4-0	2-0
Advanced Side Cutting Spatula Diameter - 0.28mm Length - 8.0mm 3/8 Circle - 135° 	<b>aurosilk</b> Silk Black Braided Silicone Coated	76cm (30") 				B6105*			
		76cm (30") 				B6106*			
Conventional Cutting Diameter - 0.43mm Length - 15.9mm 3/8 Circle - 135° 	<b>aurosilk</b> Silk Black Braided Silicone Coated	76cm (30") 						B4401*	
Taper Point Diameter - 0.66mm Length - 26.2mm 1/2 Circle - 170° 	<b>aurobond</b> Polyester Green Braided	90cm (35") 							B2202*
Taper Point Diameter - 0.71mm Length - 17.5mm 1/2 Circle - 180° 	<b>aurobond</b> Polyester Green Braided	90cm (35") 							E2204*
Reverse Cutting Diameter - 0.38mm Length - 13.1mm 3/8 Circle - 130° 	<b>aurolene</b> Polypropylene Blue Monofilament	70cm (28") 					P5301*		
Taper Point Diameter - 0.10mm Length - 5.1mm 3/8 Circle - 135° 	<b>aurolene</b> Polypropylene Blue Monofilament	15cm (6") 	P9201*						
Advanced side cutting Spatula Diameter - 0.33 mm Length - 8.0 mm 1/4 Circle - 90° 	<b>Polycryl</b> PolyGlycolic Acid Violet Braided	45cm (18") 				A6108			
Advanced Micro Point Spatula Diameter - 0.20 mm Length - 6.0 mm 3/8 Circle - 140° 	<b>Polycryl</b> PolyGlycolic Acid Violet Braided	30cm (12") 		A8104					
Advanced Micro Point Spatula Diameter - 0.230mm Length - 6.5 mm 3/8 Circle - 140° 	<b>Polycryl</b> PolyGlycolic Acid Violet Braided	30cm (12") 			A7102				
Taper Point Diameter - 0.430mm Length - 16mm 1/2 Circle - 180° 	<b>aurosilk</b> Silk Black Braided Silicone Coated	90cm (36") 						B4203	

Information published in this catalog is subject to change without notification